PENGUKUR VOLUME AKTUAL TANAH.

Judul Asli : SOIL ACTUAL VOLUMENOMETER.

Alih bahasa

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Sulusuban, April 1993.

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331

331 SOIL ACTUAL VOLUMENOMETER  
Model 100

This instrument measures the actual volume (volumetric sum of solid phase and liquid phase) of the soil sample collected or filled up by sampling cylinder with whole volume of 100 ml.

Sample capacity: 100 ml  
Measuring range: 25 - 100 ml  
Accuracy:  $\pm 0.1$  ml  
Dimensions: 40 x 17 x 55 cm  
Weight: 12 kg  
Accessories:  
Sample cylinders .... 1 set  
(6 pcs. with case)

INSTRUCTIONS

---

№ 331-B

55

KIYA SEISAKUSHO, LTD.

---

# Kiya Seisakusho, Ltd.

Manufactures and Distributors of Agricultural and General Laboratory Equipment

20-8, Mukogaoka 1-chome, Bunkyo-ku, Tokyo 113, Japan.

Phone: (03) 813-1221

Kiya No. 331-B




## SOIL ACTUAL VOLUMENOMETER SET

Model 100S

This is a set consists of necessary tools for soil survey by using the actual volumenometer(Model 100).

### \* Components and Specifications

- 1) Soil actual volumenometer .....1 unit
  - Sample capacity: 100 ml
  - Measuring range: 25 - 100 ml
  - Accuracy: within 0.1 ml
  - Dimensions: 40 x 17 x 55 cm
  - Weight: 12 kg
- 2) Soil sample cylinders .....10 sets
  - Capacity: 100 ml
  - Materials: made of nickel plated brass
  - No. of a set: 6 pcs.
  - Case: canvas, hand carrying type
- 3) Soil sampler .....1 set
  - Sampling capacity: 100 ml
  - Sampling depth: 100 cm
  - (with extention rods, 2 pcs. of 50 cm)
  - Case: canvas, shoulder type
  - Accessories: brush and knife
- 4) Soil auger .....1 set
  - Edge diameter: 10 cm
  - Auger depth: 100 cm
  - (with extension rods, 2 pcs. of 50 cm)
  - Case: canvas, shoulder type
  - Accessories: 2 spanners

# Kiya Seisakusho, Ltd.

Manufactures and Distributors of Agricultural and General Laboratory Equipment

20-8, Mukogaoka 1-chome, Bunkyo-ku, Tokyo 113, Japan.

Phone: (03) 513-1221

Kiya No. 331

## SOIL ACTUAL VOLUMENOMETER


Model 100

(Three Phase Distribution Measuring Instrument)

This instrument measures the actual volume (volumetric sum of solid phase and liquid phase) of the soil sample collected or filled up by sampling cylinder with whole volume of 100 ml. Three phase distribution of soils such as moisture ratio, porosity, and solid phase ratio can be derived from the relation between the said value (V) the weight (W) and the specific gravity (d). It is noteworthy that the measurement can be completed within 1 or 2 minutes with high accuracy without any thermal dehydration.

### \* Specifications

Sample capacity:	100 ml
Measuring range:	25 - 100 ml
Accuracy:	$\pm 0.1$ ml
Dimensions:	40 x 17 x 55 cm
Weight:	12 kg
Accessories:	Sample cylinders (6 pcs. with case) .....1 set

331

• Principle

- 1) The soil is a substance having a grained structure, and consists of three phases, i.e., the solid phase, the liquid phase and the air phase, as shown in Fig.1. The actual volume means the total of the volumes of the solid phase and the liquid phase. This unit is an instrument measuring the actual volume of the soil directly.
- 2) When the actual volume is known, the fundamental physical properties of the soil can be obtained very simply and systematically one by one, as shown in Fig.3, such a method is called the actual volume method.
- 3) The instrument has a structure shown in Fig.2; measuring tubes and U-shaped tubes are all made of acrylic resin, and the metal parts are of alloy or brass, and they are correlated. The measuring range is from 25 ml to 100 ml, and the accuracy is  $\pm 0.1$  ml.
- 4) The sample cylinder can be commonly used for not only the actual volumeter but a soil pF measuring apparatus and a soil permeability measuring apparatus; it has an inner diameter of 50 mm, a height of 51 mm, and the inner capacity of 100 ml; it is made of metal and has the weight of about 80 g. For collecting the soil in nature, a soil sampler(optional parts) is used. (Refer to Fig.4)

• Operation

- 1) Turn the cock K to the left to open.
- 2) Open the cock K1 to adjust the liquid surface at the center of the measuring tube B, and close K1.
- 3) Open the cock K3 to adjust the liquid surface at the point of 60 ml in the measuring tube A, and close K3.
- 4) Open the cock K2 to adjust the liquid surface at the measuring point in the lower part of the U-shaped tube, and close K2.


Fig. 1


Fig. 2

66P

- 5) Remove the cover C of the sample chamber, put a basin, a net plate, a sample cylinder, and test pieces of 60 ml (10, 20, 30 ml) of the instrument into, and close the chamber tightly.
- 6) Turn the cock K to the right to close, open the cock K2 to raise the liquid surface in the U-shaped tube, and when the surface reaches the upper point E, close K2.
- 7) Adjustment of point O is completed when all of the liquid surfaces are at point E.
  - a) When the right liquid surface is higher, lower the liquid surface of the micro tube B.
  - b) When the left liquid surface is higher, raise the liquid surface of the micro tube B.

\* Note: This operation is performed by opening K2, lowering the liquid surface in the U-shaped tube to the lower observation point 1, and opening K.
- 8) After adjusting point O, adjust point O of the sliding scale L to the liquid surface of the micro tube B.
- 9) Open the cock K2, adjust the liquid surface of the U-shaped tube to the lower observation point 1, close K2, and open the cock K.
- 10) Take the test pieces and the sample cylinder out from the sample chamber and set the collected soil as removing the cover of the cylinder. At this time, confirm whether the rubber packing of the sample chamber is clean or not.
- 11) Open the cock K3, set the liquid surface of the measuring tube A at the predicted value, and close K3. Close the cock K, open the cock K2, and raise the liquid surface up to the upper point E in the U-shaped tube. When the difference of the liquid surface at the upper point is large, adjust the liquid surface of the measuring A and when it is small, adjust the liquid surface of the tube B.

\* Note: When lowering or raising the liquid surface of the measuring tube A or B, be sure to lower the liquid surface of the U-shaped tube, open K and then adjust.
- 12) When the liquid surface is fit at the upper point E in the U-shaped tube, the value is the actual volume of the soil in the sample chamber.

$$V = A + B \quad \text{or} \quad V = A - B$$

where; V: actual volume of soil sample

A: measuring tube A

B: measuring tube B

\* Note 1: Don't open the cock K, the sample chamber or cocks of K1 and K3 when keeping the liquid surface of the U-shaped tube raised (the whole system is under a compressed state).

\* Note 2: When measuring, avoid the direct sunshine as much as possible (for keeping the whole system at an equal temperature as much as possible).

• Soil Structure and Term

1) Total volume	$V_t$
2) Total weight	$W$
3) Actual volume	$V$
4) Weight of solid phase	$S = W - M$
5) Volume of solid phase	$V_s = (W - V)/(d - 1)$
6) Weight of soil water	$M = V_l$
7) Volume of soil water	$V_l = V - V_s$
8) Volume of soil air	$V_a = 100 - V$
9) Porosity	$P = 100 - V_s$
10) Solid ratio	$S_v = V_s$
11) Water ratio	$M_v = V_l$
12) Air ratio	$A = V_a$
13) Water saturation percentage	$H = M_v/P \times 100$
14) Air percentage	$U = 100 - H$


Fig. 3

15) Water content by wet soil	$M_m = M/W \times 100$
16) Water content by weight	$M_o = M/S \times 100$
17) Specific gravity	$d = S/V_s$
18) Apparent specific gravity or Bulk density	$d_o = S/100$
19) Actual specific gravity	$d_m = W/V$
20) Water-solid ratio	$L_s = V_l/V_s$


Fig. 4: Soil Sampler & Cylinders